## INSTRUCTIONS FOR USE

EN

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FR

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IT

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PL

1. Remove joint effusion, if present.
2. Twenty to thirty minutes before use, remove the product box from the refrigerator, remove the blister pack from the box and allow the syringe to come to room temperature.
3. **Peel off the blister backing.** (Fig. 1)
4. While holding the blister open side down, bend the blister and allow the syringe to fall gently onto a clean surface. (Fig. 2)  
Alternatively, hold the blister open side up and bend back the blister until the barrel's luer end is exposed. Gripping the luer end of the barrel, remove the syringe from the blister.  
**Do not remove the syringe from the plunger end.**
5. **Unscrew the tip cap from the syringe** (Fig. 3.1, 3.2) **and attach an appropriately sized sterile needle**, for example 17 to 21 gauge. (Fig. 4.1, 4.2)  
**Attention:** Do not apply pressure to the plunger rod while the needle is being affixed. Verify that the needle is properly locked to the Luer Lock Adaptor (LLA). Do not overtighten the LLA; this can lead to loosening of the LLA from the barrel.
6. Apply gentle pressure to the plunger in order to expel air from the syringe needle and to verify that the syringe is operating properly. (Fig. 5)
7. The syringe is ready for use.
8. Inject intra-articularly into the joint synovial capsule using strict aseptic injection procedures.
9. For single use only. Do not re-use or resterilize.
10. Store refrigerated at 2°-8°C. Do not freeze. May be removed from refrigeration within the product shelf life and stored at room temperature (15°-25°C) for up to six (6) months.
11. Three treatments, each of 2 ml per knee should be administered, each one week apart.

Safe disposal of device: The syringe and needle should be disposed of as is done for other similar medical products. Used syringes may be disposed of with regular trash or as per glass waste. Needles may be disposed of with regular trash or as per sharps waste.

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
Fig. 4.2

Fig. 5

|                                                                                     |                                                                        |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------|
|  | Manufacturer                                                           |
|  | Authorized representative in the European Community/<br>European Union |
|  | Importer                                                               |
|  | Date of manufacture                                                    |
|  | Use-by date                                                            |
|  | Batch code                                                             |
|  | Sterilized using aseptic processing techniques                         |
|  | Do not re-sterilize                                                    |
|  | Do not use if package is damaged and consult instructions<br>for use   |
|  | Single sterile barrier system with protective packaging<br>outside     |
|  | Keep away from sunlight                                                |
|  | Keep dry                                                               |
|  | Temperature limit                                                      |
|  | Do not re-use                                                          |
|  | Consult instructions for use                                           |
|  | Caution                                                                |
|  | Medical device                                                         |
|  | Unique device identifier                                               |
|  | Complies with the provisions of the applicable EU legislation          |

 Bio-Technology General (Israel) Ltd.  
Be'er Tuvia Industrial Zone, POB 571, Kiryat Malachi, 8310402,  
Israel  
Tel: +972-8861-2020

EU REP

Ferring Pharmaceuticals A/S International PharmaScience Center (IPC)  
Amager Strandvej 405, 2770 Kastrup, Denmark  
Tel: 00 45 88 33 88 34  
Fax: 00 45 88 33 88 80

EUFLEXXA is a trademark of Ferring BV

A copy of the summary of safety and clinical performance (SSCP) is available upon request from the manufacturer. It will be found in the EUDAMED database (<https://ec.europa.eu/tools/eudamed>) when the database is fully functional.

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